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Cleaning / maintenance instruction Ear-syringe 50cc –100cc and 150cc.

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IMPORTANT

We recommend to perform undermentioned cleaning / maintenance after each use because wastage can damage the instrument and harden the cleaning process. Consistent maintenance will also considerably lengthen the instruments durability. The cleaning will take only a few minutes.

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After cleaning, each of the components has to be abundantly rinsed, this is important to make sure that no soap residues remain.

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- **Apply a few drops of new lubrication on the silicone O-ring (see photo below).**
Good lubrication is essential for the maintenance and functioning of the ear-syringe.

We recommend using a synthetic, bio-degradable lubricant specifically designed for application on medical instruments. The advantage compared to traditional oils and grease is the absence of toxic substances (therefore it is not harmful to humans), it is not degraded by contact with water, and it provides optimal lubrication without harming any materials (such as the silicone O-ring).


IMPORTANT:

To guarantee correct functioning of the ear syringe, it is essential that this step is performed. Proper lubrication makes sure that the piston can move through the tube freely (without too much friction). If the syringe is not properly lubricated, the O-ring might get off the Piston during treatment (because of too much friction). This creates a risk for the patient.

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- **Reassemble the ear-syringe.**
Place back the piston slightly slant, to protect the silicone O-ring (see photo below). Then gently (to avoid damaging the threads) twist the end cap on to the tube. Finally, move the piston up and down several times to evenly spread the lubrication and to check the functioning of the ear-syringe.

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The recommended measures should be implemented in accordance with the applicable hygiene rules and guidelines for workplace safety.