



Vscan AirTM

See More. Treat Faster.



Vscan



Ready for a new view?

Get crystal-clear images with the power of a high-performance ultrasound machine in a small, lightweight and wireless solution.

Insight right in your pocket

Supports early diagnosis and helps you quickly decide on next steps.

Developed for clinicians

It is easy to use, intuitive to navigate and optimized for your mobile device**



Dual-probe imaging power

Complete both shallow and deep scans without switching probes or compromising on image quality.

Durability, defined

Withstands everyday handling from scratches to drops and even a rainstorm.*



*Probe withstands drops: tested to Military drop standards (MIL-810G).
1P67 waterproof, compatible with high level disinfection techniques.

**Works with a range of Android and iOS smartphone and tablets.

Experience handheld ultrasound without compromising on portability or image quality



Break free.

Take Vscan Air anywhere* with wireless freedom and maximum portability. Its small, lightweight and wireless dual-probe provides you with a convenient, durable and flexible handheld ultrasound solution.



Clear images.
Informed decisions.


Engineered with SignalMax™ technology, Vscan Air delivers sharp images to support early disease detection and help monitor treatment.



By your side.

We support your ultrasound journey from first scan to confident use. Make the most of your system from the easy-to-use probe and intuitive app to GE Healthcare's renowned service and educational tools. Vscan Air comes with a 3-year worry-free warranty.

*The device has been verified for limited use outside of professional healthcare facilities. Use is restricted to environmental properties described in the user manual. Please contact your GE Healthcare sales representative for detailed information.



Securely store, share and collaborate

Flexible wireless dataflow:

- DICOM®, including query from Modality Worklist server, Store to DICOM PACS, receive Storage Commitment
- JPG/MPG export to shared network folders (supports importing into EMR)
- Share auto-anonymized images with other apps
- No cloud required



DICOM server



Shared folder



Other apps

Security:

- Secure DICOM including the option of Peer Authentication
- Secure User Authentication
- Passwords required to access patient data

SignalMax™ Technology

Power of a high-performance system
in a pocket-sized solution

See into the whole body with
dual-probe (2 in 1) technology

Console-quality signal processing
in a miniaturized chipset

Stunning image quality in a
wireless handheld device

Builds upon GE's expertise in
ultrasound innovation

Vscan Air™
See More. Treat Faster.



© 2021 General Electric Company – All rights reserved. GE, the GE Monogram and Vscan Air are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc. doing business as GE Healthcare.

DICOM is a registered trademark of the National Electrical Manufacturers Association.

Products mentioned in the material may be subject to government regulation and may not be available in all countries. Shipment and the effective sale can only occur after approval from the regulator. Please check with your local GE Healthcare representative for details.

JB02242XX